



## WORLD MEDICO | Sensitive

- Blå Sensitive Nitril
- Godkendt til kontakt med fødevarematerialer
- Godkendt til medicinsk brug
- Puddorfri



KATEGORI	EN-STANDARD
 CAT III	EN   

MATERIALER	FORPAKNING
100% Nitril	100 stk. pr. æske

Vare nr.	EAN	BRAND	MODEL	STØRRELSER
10007	6975390730728	WORLD	SENSITIVE	Small
10008	6975390730735	WORLD	SENSITIVE	Medium
10009	6975390730742	WORLD	SENSITIVE	Large
10010	6975390730759	WORLD	SENSITIVE	XL

Disse første classes 4 grams handsker, giver en unik beskyttelse mod bla. kemikalier. Handskerne sikrer høj komfort og pasform. Undersøgelseshandskerne meget slidstærke men er samtid ekstremt bløde og elastiske . De tilpasser sig håndens form under brug, og giver en god gennemfølgelse. Handskerne er åndbare men uigennemtrængelige. Handskerne sikrer et solidt greb og kan bruges til en lang række opgaver såsom rengøringsopgaver, plejeopgaver, håndtering af fødevarer, laboratorieopgaver, medicinhåndtering, cytotatika og hormoncremer. Indeholder ingen puder eller latexproteiner og kan derfor sikkert anvendes ved latexallergi

God og kraftig engangshandske med rullekant.

Tykkelse:

Håndflade: 0,07 mm (+/- 0,02)

Fingre: 0,09 mm (+/- 0,02)

Manchet: 0,05 mm (+/- 0,02)

### Leverandør

Firma: World Medico ApS

Adresse: Østerfælled Torv 3

Post nr. / By: 2100 / København Ø

Land: Danmark

Telefon / E-mail: +45 53692712 / info@worldmedico.dk

<b>CE2777</b> <b>CAT III</b>	These products are complied as Medical Device Class I, Category III Personal Protective Equipment (PPE) by the European PPE Regulation 2016/425 and have been shown to comply with this Regulation through the Harmonised European Standard(s): EN ISO 21420:2020, EN ISO 374-1:2016+A1:2018, EN ISO 374-5:2016		<b>EN ISO 374-1:2016/Type B</b>  <b>KFT</b> <b>EN ISO 374-5:2016</b>  <b>VIRUS</b> Resistance to viruses, bacteria and Fung																					
	<b>EN ISO 374-1:2016+A1:2018 Permeation Levels are based on breakthrough times as follows:</b> <table border="1"> <thead> <tr> <th>Chemical</th> <th>Permeation Level</th> </tr> </thead> <tbody> <tr> <td>&gt; 10min</td> <td>1</td> </tr> <tr> <td>&gt; 20min</td> <td>2</td> </tr> <tr> <td>&gt; 40min</td> <td>3</td> </tr> <tr> <td>&gt; 120min</td> <td>4</td> </tr> <tr> <td>&gt; 240min</td> <td>5</td> </tr> <tr> <td>&gt; 480min</td> <td>6</td> </tr> </tbody> </table>	Chemical	Permeation Level	> 10min	1	> 20min	2	> 40min	3	> 120min	4	> 240min	5	> 480min	6	<b>ISO 374-5:2016, class 1.</b> The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.  Before usage, inspect the gloves for any defector imperfections. Components used in glove manufacturing may cause allergic reactions. If reactions occur then seek medical advice immediately. Only use with dry and clean hands.	<b>EN ISO 374-1:2016+A1:2018 /Type B</b> <table border="1"> <thead> <tr> <th>Level</th> <th>EN ISO 374-4:2019 Degradation (%)</th> </tr> </thead> <tbody> <tr> <td>6</td> <td>-5,9</td> </tr> <tr> <td>2</td> <td>-26,4</td> </tr> <tr> <td>4</td> <td>-2,1</td> </tr> </tbody> </table>	Level	EN ISO 374-4:2019 Degradation (%)	6	-5,9	2	-26,4	4
Chemical	Permeation Level																							
> 10min	1																							
> 20min	2																							
> 40min	3																							
> 120min	4																							
> 240min	5																							
> 480min	6																							
Level	EN ISO 374-4:2019 Degradation (%)																							
6	-5,9																							
2	-26,4																							
4	-2,1																							
<b>EU Type Examination and On-going Conformity Module C2, Carried out by Noti =d Body 2777, SAKKA Technology Europe Ltd, Bracetown Business Park, Clonsilla, Dublin 15, D15 3K2P Ireland</b>	<b>Declaration of Conformity</b>		<b>EN ISO 374-4:2019</b> Degradation results indicate the change in puncture resistance of the gloves after exposure to the challenge chemical.																					




**worldmedico**  
 We Supply with Passion  
**NITRILE | EXAMINATION**  
 Gloves Blue  
 Sensitive  
 www.worldmedico.eu

**10 Boxes x 100 pcs**

**MODEL: 10-STRT01**  
**STORAGE METHOD:**  
 The package shall be stored in room with normal temperature lower than 40°C, please keep away from the sun, and avoid direct contact with fire.  
**SPECIFICATIONS:** 10 Examination  
**EXECUTIVE STANDARD:** EN 455  
**PACKAGING SPECIFICATIONS:** 10\*10\*10cm  
**SHELF LIFE:** 3 years

**S**


**worldmedico**  
 We Supply with Passion  
**NITRILE | EXAMINATION**  
 Gloves Blue  
 www.worldmedico.eu



**Leverandør**  
 Firma: World Medico ApS  
 Adresse: Østerfælled Torv 3  
 Post nr. / By: 2100 / København Ø  
 Land: Danmark  
 Telefon / E-mail: +45 53692712 / info@worldmedico.dk

